

# Product Datasheet

## Anti-Human CD22 Monoclonal Antibody, Unconjugated Clone MEM-01 853.193.020

Product Name	Anti-Human CD22 Monoclonal Antibody, Unconjugated Clone MEM-01
Catalog Number	853.193.020
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Purified
Tested Applications	FC, IP, WB
Immunogen	Burkitt's lymphoma cell line Raji
Brand	Diaclone
Form/Appearance	Phosphate buffered saline (PBS) with 15 mM sodium azide, approx. pH 7.4
Storage	Store at +2-8°C. DO NOT FREEZE.
Note	We recommend that each laboratory determine an optimum working titre for use in its particular application
Isotype	IgG1
Clonality	Monoclonal
Specificity	Reacts with CD22 (BL-CAM), a 130 kDa type I transmembrane glycoprotein (immunoglobulin superfamily) expressed in the cytoplasm of pro-B and pre-B lymphocytes