

Product Datasheet

Ferritin, Human Liver 270-40

Product Name	Ferritin, Human Liver
Catalog Number	270-40
Alternative Names	FTL
Brand	Lee Biosolutions
Form/Appearance	Solution in tris buffered sodium chloride, pH 7.5, with bromo-nitro-dioxane and methylisothiazolone as preservatives. Clear, reddish-brown solution
Concentration	3 - 6 mg/mL (Lowry)
Storage	+2-8 °C
Note	Contains both heavy chains (H) and light chains (L)
Purity	≥ 95% (SDS-PAGE)
CAS Number	9007-73-2
Source	Human Liver
Molecular Weight	440000
Shelf Life	3 years