

Product Datasheet

Ferritin, Human Spleen 270-50

Product Name	Ferritin, Human Spleen
Catalog Number	270-50
Brand	Lee Biosolutions
Form/Appearance	Liquid Clear, reddish-brown solution
Concentration	2.0 - 7.0 mg/mL (Lowry)
Storage	+2-8 °C
Note	Contains both heavy chains (H) and light chains (L)
Purity	≥ 95% (SDS-PAGE)
CAS Number	9007-73-2
Source	Human Spleen
Molecular Weight	440000
Shelf Life	2 years