Medix Biochemica

Oy Medix Biochemica Ab

Headquarters Klovinpellontie 3 Fl-02180 Espoo

Finland

Email: sales@medixbiochemica.com

Medix Biochemica USA Inc.

10850 Metro Court Maryland Heights, MO 63043 United States of America

Email: sales@medixbiochemica.com

Product Datasheet

Anti-Human CD30 Monoclonal Antibody, Unconjugated Clone MEM-268 853.203.020

Product Anti-Human CD30 Monoclonal Antibody, Unconjugated Clone

Name MEM-268

Catalog 853.203.020 Number

Description Mouse monoclonal

Species/Host Mouse

Conjugation Purified

Tested IHC, FC **Applications**

Immunogen Expression vector containing CD30 cDNA (booster suspension of

THP-1 cell line)

Brand Diaclone

Form/AppearaPhosphate buffered saline (PBS) with 15 mM sodium azide,

approx. pH 7.4

Storage Store at +2-8°C. DO NOT FREEZE.

Note We recommend that each laboratory determine an optimum

working titre for use in its particular application

Isotype IgG

Clonality Monoclonal

Specificity Recognizes extracellular part of CD30 (Ki-1 antigen), a 105 kDa

single chain glycoprotein expressed on Hodgkin's and

Reed-Sternberg cells