

Product Datasheet

Anti-Human CD30 Monoclonal Antibody, Unconjugated Clone MEM-268 853.203.020

Product Name	Anti-Human CD30 Monoclonal Antibody, Unconjugated Clone MEM-268
Catalog Number	853.203.020
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Purified
Tested Applications	FC, IHC
Immunogen	Expression vector containing CD30 cDNA (booster suspension of THP-1 cell line)
Brand	Diaclone
Form/Appearance	Phosphate buffered saline (PBS) with 15 mM sodium azide, approx. pH 7.4
Storage	Store at +2-8°C. DO NOT FREEZE.
Note	We recommend that each laboratory determine an optimum working titre for use in its particular application
Isotype	IgG
Clonality	Monoclonal
Specificity	Recognizes extracellular part of CD30 (Ki-1 antigen), a 105 kDa single chain glycoprotein expressed on Hodgkin's and Reed-Sternberg cells