

Product Datasheet

Haptoglobin, Porcine Serum 333-11

Product Name	Haptoglobin, Porcine Serum
Catalog Number	333-11
Alternative Names	Hp
Brand	Lee Biosolutions
Form/Appearance	Frozen solution
Concentration	≥ 0.5 mg/mL (A ₂₈₀)
Storage	−20 °C
Purity	≥ 90% (SDS-PAGE)
Source	Porcine Serum
Buffer	Phosphate buffered saline with 0.05% sodium azide
Shelf Life	2 years