

Product Datasheet

Haptoglobin, Mouse Serum 334-11

Product Name	Haptoglobin, Mouse Serum
Catalog Number	334-11
Alternative Names	Hp
Brand	Lee Biosolutions
Form/Appearance	Liquid Frozen solution
Concentration	≥ 0.5 mg/mL (Coomassie)
Storage	-20 °C
Purity	≥ 90% (SDS-PAGE)
CAS Number	9087-69-8
Source	Mouse Serum
Molecular Weight	alpha Hp: ~10,000 beta Hp: ~40,000
Shelf Life	2 years