

Product Datasheet

Haptoglobin, Bovine Serum 337-11

Product Name	Haptoglobin, Bovine Serum
Catalog Number	337-11
Alternative Names	HAP, HP
Brand	Lee Biosolutions
Form/Appearance	Liquid
Concentration	≥ 0.5 mg/mL (A ₂₈₀)
Storage	-20 °C
Note	Due to its high affinity for hemoglobin, haptoglobin preparations unavoidably contain hemoglobin, which co-purifies during all chromatographic steps
Purity	$\geq 90\%$ (SDS-PAGE)
CAS Number	9087-69-8
Source	Bovine Serum
Buffer	Phosphate buffered saline with 0.05% sodium azide
Shelf Life	2 years