

Product Datasheet

Hemopexin, Human Plasma - High Purity 339-18

Product Name	Hemopexin, Human Plasma - High Purity
Catalog Number	339-18
Alternative Names	Haemopexin, HPX, beta-1B-glycoprotein
Brand	Lee Biosolutions
Form/Appearance	Lyophilized
Concentration	Reported in mg/vial ($E^{0.1\%}_{280nm} = 1.97$)
Storage	-20 °C
Purity	≥ 95% (SDS-PAGE)
CAS Number	9013-71-2
Source	Human Plasma
Buffer	Lyophilized from phosphate buffered saline, pH 7.4
Molecular Weight	~57,000
Shelf Life	5 years