

## Product Datasheet

### Immunoglobulin A (IgA), Human Plasma 340-12

<b>Product Name</b>	Immunoglobulin A (IgA), Human Plasma
<b>Catalog Number</b>	340-12
<b>Alternative Names</b>	IgA
<b>Brand</b>	Lee Biosolutions
<b>Form/Appearance</b>	Lyophilized White to off-white powder
<b>Concentration</b>	Reported in mg/vial ( $A_{280}$ , $E^{1\text{mg/mL}} = 1.34$ )
<b>Storage</b>	-20 °C
<b>Purity</b>	≥ 96% (SDS-PAGE)
<b>Source</b>	Human Plasma
<b>Shelf Life</b>	5 years