

Product Datasheet

Immunoglobulin A (IgA), Human Plasma 340-13

| | |
|--------------------------|---------------------------------------------------|
| Product Name | Immunoglobulin A (IgA), Human Plasma |
| Catalog Number | 340-13 |
| Alternative Names | IgA, sIgA |
| Brand | Lee Biosolutions |
| Form/Appearance | Liquid Clear, colorless to slight yellow solution |
| Concentration | 3 - 6 mg/mL (A_{280} , $E^{0.1\%} = 1.32$) |
| Storage | 2-8° C |
| Purity | ≥ 95% (SDS-PAGE) |
| Source | Human Plasma |
| Specificity | Mucosal Immunity |
| Molecular Weight | 160000 |
| Shelf Life | 2 years |