

## Product Datasheet

### Immunoglobulin G Fab (IgG Fab) Antigen from human plasma 340-27

<b>Product Name</b>	Immunoglobulin G Fab (IgG Fab) Antigen from human plasma
<b>Catalog Number</b>	340-27
<b>Alternative Names</b>	IgG Fab, Fragment antigen-binding
<b>Brand</b>	Lee Biosolutions
<b>Form/Appearance</b>	Liquid Clear, colorless liquid
<b>Concentration</b>	$\geq 5.0$ mg/mL ( $E^{1\text{ mg/mL}}_{280\text{nm}} = 1.35$ )
<b>Storage</b>	-20 °C
<b>Purity</b>	$\geq 95\%$ (SDS-PAGE)
<b>Source</b>	Human Plasma
<b>Molecular Weight</b>	~50,000
<b>Shelf Life</b>	3 years