

Product Datasheet

Immunoglobulin G Fab (IgG Fab) Antigen from human plasma 340-27

Product Name	Immunoglobulin G Fab (IgG Fab) Antigen from human plasma
Catalog Number	340-27
Alternative Names	IgG Fab, Fragment antigen-binding
Brand	Lee Biosolutions
Form/Appearance	Liquid Clear, colorless liquid
Concentration	≥ 5.0 mg/mL ($E^{1\text{ mg/mL}}_{280\text{nm}} = 1.35$)
Storage	-20 °C
Purity	$\geq 95\%$ (SDS-PAGE)
Source	Human Plasma
Molecular Weight	~50,000
Shelf Life	3 years