

Product Datasheet

Immunoglobulin M (IgM) from human myeloma plasma 340-35

Product Name	Immunoglobulin M (IgM) from human myeloma plasma
Catalog Number	340-35
Alternative Names	IgM
Brand	Lee Biosolutions
Form/Appearance	Liquid Clear, colorless liquid
Concentration	1.0 - 5.0 mg/mL ($E^{0.1\%}_{280nm} = 1.18$)
Storage	-20 °C
Purity	≥ 95% (SDS-PAGE)
Source	Human Myeloma Plasma
Molecular Weight	950000
Shelf Life	2 years