

# Product Datasheet

## Anti-Human CD41a Monoclonal Antibody, Unconjugated Clone HIP8 853.223.020

Product Name	Anti-Human CD41a Monoclonal Antibody, Unconjugated Clone HIP8
Catalog Number	853.223.020
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Purified
Tested Applications	IHC, FC
Brand	Diaclone
Form/Appearance	Phosphate buffered saline (PBS) with 15 mM sodium azide, approx. pH 7.4
Storage	Store at +2-8°C. DO NOT FREEZE.
Note	We recommend that each laboratory determine an optimum working titre for use in its particular application
Isotype	IgG1
Clonality	Monoclonal
Specificity	Reacts with CD41a antigen. CD41 (GPIIb) is a transmembrane glycoprotein (integrin family) composed of two chains GPIIbα and GPIIbβ