

# Product Datasheet

## Immunoglobulin E (IgE) - Human Myeloma Plasma 340-40

|                     |  |
|---------------------|--|
| Product Name        | Immunoglobulin E (IgE) - Human Myeloma Plasma  |
| Catalog Number      | 340-40   |
| Alternative Names   | IgE  |
| Brand               | Lee Biosolutions   |
| Form/Appearance     | Liquid   |
| Storage             | -20 °C   |
| Source              | Human Myeloma Plasma   |
| Shelf Life          | Due to differing stabilities of various biomarkers in this complex biological material it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time. |
| Biological Activity | ≥ 1,500,000 U/mL (ELISA)   |