

Product Datasheet

Epididymis Protein 4 (HE4) antigen 345-50

Product Name	Epididymis Protein 4 (HE4) antigen
Catalog Number	345-50
Species/Host	E. coli
Alternative Names	HE4, WAP four-disulfide core domain protein, WFDC2
Brand	Lee Biosolutions
Form/Appearance	Liquid
Concentration	≥ 0.5 mg/mL (Immunoassay)
Storage	-20 °C
Purity	≥ 85% (SDS-PAGE)
Source	Recombinant
Molecular Weight	10036
Shelf Life	18 months
Biological Activity	Immunologically active
Tag	N-terminal His Tag