

## Product Datasheet

### 0.2 micron Filtered LDL Low Density Lipoprotein, Human Plasma 360-12

<b>Product Name</b>	0.2 micron Filtered LDL Low Density Lipoprotein, Human Plasma
<b>Catalog Number</b>	360-12
<b>Alternative Names</b>	LDL, LDL-C, Beta-Lipoprotein, Concentrate, Fraction
<b>Brand</b>	Lee Biosolutions
<b>Form/Appearance</b>	Liquid Yellow to orange liquid
<b>Storage</b>	-60 °C or below
<b>Note</b>	EDTA-free
<b>Purity</b>	Plasma Fraction
<b>Source</b>	Human Plasma
<b>Buffer</b>	Sodium chloride and sucrose
<b>Shelf Life</b>	2 years
<b>Content</b>	3,000 - 4,000 mg/dL (Enzymatic) <b>Inquire for custom concentrations and cholesterol ratios</b>
<b>Triglyceride Ratio</b>	$\leq 0.8$ (Triglycerides/Total Cholesterol)