

# Product Datasheet

## Anti-Human CD43 Monoclonal Antibody, Unconjugated Clone MEM-59 853.243.020

Product Name	Anti-Human CD43 Monoclonal Antibody, Unconjugated Clone MEM-59
Catalog Number	853.243.020
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Purified
Tested Applications	FC, IP, WB, IHC
Immunogen	Human T lymphocytes
Brand	Diaclone
Form/Appearance	Phosphate buffered saline (PBS) with 15 mM sodium azide, approx. pH 7.4
Storage	Store at +2-8°C. DO NOT FREEZE.
Note	We recommend that each laboratory determine an optimum working titre for use in its particular application
Isotype	IgG1
Clonality	Monoclonal
Specificity	Recognizes neuraminidase-sensitive epitope on CD43 (Leukosialin), a 95-135 kDa type I transmembrane glycoprotein (mucin-type) which is involved in lymphocyte activation