

Oy Medix Biochemica Ab

Headquarters Klovinpellontie 3 Fl-02180 Espoo

**Finland** 

Email: sales@medixbiochemica.com

**Medix Biochemica USA Inc.** 

10850 Metro Court Maryland Heights, MO 63043 United States of America

Email: sales@medixbiochemica.com

## **Product Datasheet**

## Anti-Human CD43 Monoclonal Antibody, Unconjugated Clone MEM-59 853.243.020

**Product** Anti-Human CD43 Monoclonal Antibody, Unconjugated Clone

Name MEM-59

Catalog Number 853.243.020

**Description** Mouse monoclonal

Species/Host Mouse

Conjugation Purified

Tested

WB, FC, IP, IHC

Applications

Immunogen Human T lymphocytes

**Brand** Diaclone

Form/AppearaPhosphate buffered saline (PBS) with 15 mM sodium azide,

approx. pH 7.4

**Storage** Store at +2-8°C. DO NOT FREEZE.

**Note** We recommend that each laboratory determine an optimum

working titre for use in its particular application

**Isotype** lgG1

**Clonality** Monoclonal

**Specificity** Recognizes neuraminidase-sensitive epitope on CD43

(Leukosialin), a 95-135 kDa type I transmembrane glycoprotein

(mucin-type) which is involved in lymphocyte activation