

Product Datasheet

Sterile Filtered HDL High Density Lipoprotein, Human Plasma 361-12

Product Name	Sterile Filtered HDL High Density Lipoprotein, Human Plasma
Catalog Number	361-12
Alternative Names	HDL, HDL-C, HDL Concentrate, HDL Fraction, alpha-Lipoprotein
Brand	Lee Biosolutions
Form/Appearance	Liquid Yellow to orange liquid
Storage	-60 °C or below
Note	EDTA-free
Purity	Plasma Fraction
Source	Human Plasma
Buffer	Sodium chloride and sucrose
Molecular Weight	175,000-360,000 (Tetz second Edition)
Shelf Life	2 years
Content	3,000 - 4,000 mg/dL (Enzymatic) Inquire for custom concentrations and cholesterol ratios
Triglyceride Ratio	≤ 0.6 (Triglycerides/Total Cholesterol)