

Product Datasheet

Anti-Human CD80 Monoclonal Antibody, FITC Conjugated Clone MEM-233 853.251.010

Product Name	Anti-Human CD80 Monoclonal Antibody, FITC Conjugated Clone MEM-233
Catalog Number	853.251.010
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	FITC
Tested Applications	FC
Immunogen	CD80 fusion protein
Brand	Diaclone
Form/Appearance	Phosphate buffered saline (PBS) containing 15 mM sodium azide and 0.2% (w/v) high-grade protease free Bovine Serum Albumin (BSA). No reconstitution is necessary
Storage	Store at +2-8°C protected from light. DO NOT FREEZE.
Note	For flow cytometry we recommend using 20 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension
Isotype	IgG1
Clonality	Monoclonal
Specificity	Reacts with CD80 (B7-1), a 60 kDa single chain type I glycoprotein (immunoglobulin supergene family) expressed on activated B and T lymphocytes and macrophages