

Product Datasheet

Lipoprotein(a) (Lp(a)), Human Plasma 400-42

Product Name	Lipoprotein(a) (Lp(a)), Human Plasma
Catalog Number	400-42
Alternative Names	Lp(a), LPA, AK38, APOA, LP
Brand	Lee Biosolutions
Form/Appearance	Liquid
Concentration	≥ 0.1 mg/mL (Immunoturbidimetric)
Storage	+2-8 °C
Purity	≥ 95% (Electrophoresis)
Source	Human Plasma
Buffer	Phosphate buffered saline with 0.3 mM EDTA and 0.05% sodium azide
Molecular Weight	Approx. 1,000,000
Shelf Life	1 year