

Product Datasheet

anti-Human Lipoprotein(a) [Lp(a)] Polyclonal Antibody 400-45

Product Name	anti-Human Lipoprotein(a) [Lp(a)] Polyclonal Antibody
Catalog Number	400-45
Description	Affinity Purified Polyclonal Antibody
Species/Host	Sheep
Immunogen	High purity human Lp(a)
Alternative Names	Lp(a), LPA, AK38, APOA, LP
Brand	Lee Biosolutions
Form/Appearance	Liquid
Concentration	0.5 - 2.0 mg/mL
Storage	-20 °C
Clonality	Polyclonal
Purity	Affinity Purified
Shelf Life	1 year