

Product Datasheet

Lipoprotein-associated Phospholipase A2 (Lp-PLA2) antigen 400-60

Product Name	Lipoprotein-associated Phospholipase A2 (Lp-PLA2) antigen
Catalog Number	400-60
Alternative Names	Lp-PLA2, platelet-activating factor acetylhydrolase, PAF acetylhydrolase, PLA2G7
Brand	Lee Biosolutions
Form/Appearance	Liquid
Concentration	≥ 0.5 mg/mL
Storage	-20 °C
Purity	≥ 90% (SDS-PAGE)
Source	Recombinant
Buffer	20 mM Tris-Cl, 500 mM NaCl, pH 8.0, 50% glycerol
Molecular Weight	45000
Shelf Life	18 months
Biological Activity	Immunologically active
Expression System	E. coli