

Product Datasheet

Anti-Human CD80 Monoclonal Antibody, Unconjugated Clone MEM-233 853.253.020

Product Name	Anti-Human CD80 Monoclonal Antibody, Unconjugated Clone MEM-233
Catalog Number	853.253.020
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Purified
Tested Applications	FC, IP
Immunogen	CD80 fusion protein
Brand	Diaclone
Form/Appearance	Phosphate buffered saline (PBS) with 15 mM sodium azide, approx. pH 7.4
Storage	Store at +2-8°C. DO NOT FREEZE.
Note	We recommend that each laboratory determine an optimum working titre for use in its particular application
Isotype	IgG1
Clonality	Monoclonal
Specificity	Reacts with CD80 (B7-1), a 60 kDa single chain type I glycoprotein (immunoglobulin supergene family) expressed on activated B and t lymphocytes and macrophages