Medix Biochemica

Oy Medix Biochemica Ab

Headquarters Klovinpellontie 3 Fl-02180 Espoo

Finland

Email: sales@medixbiochemica.com

Medix Biochemica USA Inc.

10850 Metro Court Maryland Heights, MO 63043 United States of America

Email: sales@medixbiochemica.com

Product Datasheet

Anti-Human CD105 Monoclonal Antibody, Unconjugated Clone MEM-226 853.263.020

Product Anti-Human CD105 Monoclonal Antibody, Unconjugated Clone

Name MEM-226

Catalog Number

853.263.020

Description Mouse monoclonal

Species/Host Mouse

Conjugation Purified

Tested Applications WB, FC, IP

Immunogen Recombinant Vaccinia virus

Brand Diaclone

Form/AppearaPhosphate buffered saline (PBS) with 15 mM sodium azide,

approx. pH 7.4

Storage Store at +2-8°C. DO NOT FREEZE.

Note We recommend that each laboratory determine an optimum

working titre for use in its particular application

Isotype IgG2a

Clonality Monoclonal

Specificity Reacts with CD105 (Endoglin), a 180 kDa type I integral

membrane homodimer glycoprotein expressed on vascular

endothelial cells