

Product Datasheet

Anti-Human CD105 Monoclonal Antibody, Unconjugated Clone MEM-226 853.263.020

Product Name	Anti-Human CD105 Monoclonal Antibody, Unconjugated Clone MEM-226
Catalog Number	853.263.020
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Purified
Tested Applications	IP, WB, FC
Immunogen	Recombinant Vaccinia virus
Brand	Diaclone
Form/Appearance	Phosphate buffered saline (PBS) with 15 mM sodium azide, approx. pH 7.4
Storage	Store at +2-8°C. DO NOT FREEZE.
Note	We recommend that each laboratory determine an optimum working titre for use in its particular application
Isotype	IgG2a
Clonality	Monoclonal
Specificity	Reacts with CD105 (Endoglin), a 180 kDa type I integral membrane homodimer glycoprotein expressed on vascular endothelial cells