

Product Datasheet

PSA-ACT Complex, Human - High Purity 498-11

Product Name	PSA-ACT Complex, Human - High Purity
Catalog Number	498-11
Alternative Names	Prostate Specific Antigen-Alpha-1-Antichymotrypsin Complex
Brand	Lee Biosolutions
Form/Appearance	Clear, colorless frozen liquid
Concentration	0.5 - 3.0 mg/mL (A ₂₈₀ , E ^{0.1%} = 1.0)
Storage	−20 °C
Purity	≥ 96% (SDS-PAGE)
Source	Human Semen & Serum
Buffer	0.2 micron filtered 0.05M sodium phosphate pH 7.5 with 0.15M sodium chloride and 0.09% sodium azide
Molecular Weight	81000
Shelf Life	5 years