

Product Datasheet

Secretory Immunoglobulin A, Human Colostrum 527-55

Product Name	Secretory Immunoglobulin A, Human Colostrum
Catalog Number	527-55
Alternative Names	slgA
Brand	Lee Biosolutions
Form/Appearance	Lyophilized
Concentration	Reported in mg/vial ($A_{278} E^{0.1\%} = 1.26$)
Storage	-20 °C
Purity	≥ 95% (SDS-PAGE)
Source	Human Colostrum
Cross Reactivity	Non-reactive with antisera to human albumin, IgD, IgE, IgG, IgM, and lactoferrin (IEP)
Molecular Weight	~385,000
Shelf Life	3 years