

Product Datasheet

Thyroglobulin (Tg), Human Thyroid - High Purity 528-11

Product Name Thyroglobulin (Tg), Human Thyroid - High Purity

Catalog Number 528-11

Alternative Names Tg

Brand Lee Biosolutions

Form/Appearance Lyophilized White to off-white powder

Concentration Reported in mg/vial (A_{280} , $E^{0.1\%}=1.04$)

Storage +2-8 °C

Purity $\geq 96\%$ (Native PAGE)

CAS Number 9010-34-8

Source Human Thyroid Gland

Buffer Lyophilized from 20 mM ammonium bicarbonate

Molecular Weight 660 000

Shelf Life 5 years