

Product Datasheet

Thyroxine Binding Globulin (TBG), Human Serum 529-13

Product Name	Thyroxine Binding Globulin (TBG), Human Serum
Catalog Number	529-13
Alternative Names	Serpin A7, T4-binding globulin, TBG
Brand	Lee Biosolutions
Form/Appearance	White lyophilized solid
Concentration	Reported in mg/vial (A_{280} , $E^{0.1\%}=0.9$)
Storage	+2-8 °C
Purity	≥ 98% (SDS-PAGE)
Source	Human Serum
Buffer	Lyophilized from 20 mM ammonium bicarbonate
Molecular Weight	55000
Shelf Life	5 years