

## Product Datasheet

### Anti-Human CD27 Monoclonal Antibody, Unconjugated Clone LT27 853.343.020

<b>Product Name</b>	Anti-Human CD27 Monoclonal Antibody, Unconjugated Clone LT27
<b>Catalog Number</b>	853.343.020
<b>Description</b>	Mouse monoclonal
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Purified
<b>Tested Applications</b>	FC
<b>Immunogen</b>	Human peripheral blood lymphocytes
<b>Brand</b>	Diaclone
<b>Form/Appearance</b>	Phosphate buffered saline (PBS) with 15 mM sodium azide, approx. pH 7.4
<b>Storage</b>	Store at +2-8°C. DO NOT FREEZE.
<b>Note</b>	We recommend that each laboratory determine an optimum working titre for use in its particular application
<b>Isotype</b>	IgG2a
<b>Clonality</b>	Monoclonal
<b>Specificity</b>	Reacts with CD27 (T14), a 50-55 kDa type I transmembrane glycoprotein (member of the TNF-receptor superfamily)