

## Product Datasheet

### Troponin T (Tn-T), Human Heart 550-21

<b>Product Name</b>	Troponin T (Tn-T), Human Heart
<b>Catalog Number</b>	550-21
<b>Alternative Names</b>	Tn-T, hTn-T, cardiac Tn-T, TnT, cTnT, cTn-T
<b>Brand</b>	Lee Biosolutions
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	$\geq 0.5$ mg/mL (A <sub>280</sub> )
<b>Storage</b>	-20 °C (up to three months), -70 °C (more than three months)
<b>Purity</b>	$\geq 98\%$ (SDS-PAGE)
<b>CAS Number</b>	77108-40-8
<b>Source</b>	Human Heart
<b>Molecular Weight</b>	~35,000
<b>Shelf Life</b>	2 years