

Product Datasheet

Troponin C (Tn-C), Human Heart
550-40

Product Name	Troponin C (Tn-C), Human Heart
Catalog Number	550-40
Alternative Names	Tn-c, hTn-c, Troponin C type 1
Brand	Lee Biosolutions
Form/Appearance	Liquid Frozen solution
Concentration	≥ 0.5 mg/mL (A ₂₈₀ , E ^{0.1%} =0.24)
Storage	-20 °C
Purity	≥ 95% (SDS-PAGE)
Source	Human Heart
Molecular Weight	~18,000
Shelf Life	5 years