

Product Datasheet

Vitamin D-Binding Protein (DBP), anti-Human Polyclonal Ab 675-12

Product Name	Vitamin D-Binding Protein (DBP), anti-Human Polyclonal Ab
Catalog Number	675-12
Species/Host	Goat
Immunogen	High purity, native human DPB
Alternative Names	anti-DBP
Brand	Lee Biosolutions
Form/Appearance	Liquid
Concentration	1.0 - 2.0 mg/mL
Isotype	IgG
Purity	Affinity Purified
Source	Polyclonal Antibody
Shelf Life	2 years