

## Product Datasheet

### Breast Milk, Pooled Donors, Lyophilized 991-01-P-L

<b>Product Name</b>	Breast Milk, Pooled Donors, Lyophilized
<b>Catalog Number</b>	991-01-P-L
<b>Description</b>	Collected at least 1 week postpartum.
<b>Tested Applications</b>	QC Materials
<b>Alternative Names</b>	Mothers milk, human milk, lactate, lactation, expressed milk, nursing milk, mother's milk, powder
<b>Brand</b>	Lee Biosolutions
<b>Form/Appearance</b>	Powder
<b>Storage</b>	2-8°C
<b>Source</b>	Pooled Human Donors
<b>Shelf Life</b>	Due to differing stabilities of various biomarkers in this complex biological material, it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time.
<b>Preservatives</b>	None added, inquire if needed