

# Product Datasheet

## Anti-Human CD41 Monoclonal Antibody, Unconjugated Clone MEM-06 853.353.020

Product Name	Anti-Human CD41 Monoclonal Antibody, Unconjugated Clone MEM-06
Catalog Number	853.353.020
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Purified
Tested Applications	FC
Immunogen	Human Leucocytes
Brand	Diaclone
Form/Appearance	Phosphate buffered saline (PBS) with 15 mM sodium azide, approx. pH 7.4
Storage	Store at +2-8°C. DO NOT FREEZE.
Note	We recommend that each laboratory determine an optimum working titre for use in its particular application
Isotype	IgG1
Clonality	Monoclonal
Specificity	Reacts with CD41 (GPIIb), a transmembrane glycoprotein (integrin family) composed of two chains GPIIbα (heavy chain; 120 kDa) and GPIIbβ (light chain; 23 kDa)