

## Product Datasheet

### Anti-Human CD41 Monoclonal Antibody, Unconjugated Clone MEM-06 853.353.020

<b>Product Name</b>	Anti-Human CD41 Monoclonal Antibody, Unconjugated Clone MEM-06
<b>Catalog Number</b>	853.353.020
<b>Description</b>	Mouse monoclonal
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Purified
<b>Tested Applications</b>	FC
<b>Immunogen</b>	Human Leucocytes
<b>Brand</b>	Diaclone
<b>Form/Appearance</b>	Phosphate buffered saline (PBS) with 15 mM sodium azide, approx. pH 7.4
<b>Storage</b>	Store at +2-8°C. DO NOT FREEZE.
<b>Note</b>	We recommend that each laboratory determine an optimum working titre for use in its particular application
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Specificity</b>	Reacts with CD41 (GPIIb), a transmembrane glycoprotein (integrin family) composed of two chains GPIIb $\alpha$ (heavy chain; 120 kDa) and GPIIb $\beta$ (light chain; 23 kDa)