

# Product Datasheet

## Pasteurized Breast Milk, Pooled Human Donors 991-01-P-PM

Product Name	Pasteurized Breast Milk, Pooled Human Donors
Catalog Number	991-01-P-PM
Alternative Names	Mothers milk, lactate, lactation
Brand	Lee Biosolutions
Form/Appearance	Liquid White to cream colored liquid
Storage	-20 °C
Source	Pooled Human Donors
Shelf Life	Due to differing stabilities of various biomarkers in this complex biological material it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time.
Preservatives	None (Available by request)