

Product Datasheet

Anti-Human CD55 Monoclonal Antibody, Unconjugated Clone MEM-118 853.373.020

Product Name	Anti-Human CD55 Monoclonal Antibody, Unconjugated Clone MEM-118
Catalog Number	853.373.020
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Purified
Tested Applications	FC, IP, WB, IHC
Immunogen	T cell line HPB-ALL
Brand	Diaclone
Form/Appearance	Phosphate buffered saline (PBS) with 15 mM sodium azide, approx. pH 7.4
Storage	Store at +2-8°C. DO NOT FREEZE.
Note	We recommend that each laboratory determine an optimum working titre for use in its particular application
Isotype	IgM
Clonality	Monoclonal
Specificity	Recognizes an epitope in SCR4 domain of CD55 (Decay accelerating factor, DAF), a 60-70 kDa glycosylphosphatidylinositol (GPI)-anchored single chain glycoprotein