

Product Datasheet

Anti-Human CD55 Monoclonal Antibody, Unconjugated Clone MEM-118 853.373.020

Product Name	Anti-Human CD55 Monoclonal Antibody, Unconjugated Clone MEM-118
Catalog Number	853.373.020
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Purified
Tested Applications	FC, IP
Brand	Diaclone
Form/Appearance	Phosphate buffered saline (PBS) with 15 mM sodium azide, approx. pH 7.4
Storage	Store at +2-8°C. DO NOT FREEZE.
Note	We recommend that each laboratory determine an optimum working titre for use in its particular application
Isotype	IgM
Clonality	Monoclonal
Source	T cell line HPB-ALL
Specificity	Recognizes an epitope in SCR4 domain of CD55 (Decay accelerating factor, DAF), a 60-70 kDa glycosylphosphatidylinositol (GPI)-anchored single chain glycoprotein
Immunisation	T cell line HPB-ALL