

Product Datasheet

Pooled Human Urine - Normal - Fresh, Never Frozen 991-03-P-FNF

Product Name	Pooled Human Urine - Normal - Fresh, Never Frozen
Catalog Number	991-03-P-FNF
Brand	Lee Biosolutions
Form/Appearance	Liquid Light to dark yellow, clear liquid
Storage	2-8° C
Purity	Raw Urine
Source	Pooled Human Donors
Shelf Life	Due to differing stabilities of various biomarkers in this complex biological material it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time.
Preservatives	None (Available by request)