

Product Datasheet

Urine, Pooled Donors, 0.2 um Filtered 991-03-P-FTD

Product Name	Urine, Pooled Donors, 0.2 um Filtered
Catalog Number	991-03-P-FTD
Tested Applications	QC Materials
Alternative Names	Urea, excretory fluid, pee, bodily waste, urinary, liquid waste, clarified, processed, filter, urinalysis
Brand	Lee Biosolutions
Form/Appearance	Liquid
Storage	-20°C
Source	Pooled Human Donors
Shelf Life	If unopened, up to 3 months at 2-8°C and >3 months at -20°C. Due to differing stabilities of various biomarkers in this complex biological material, it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time.
Preservatives	None added, inquire if needed