

Product Datasheet

Anti-Human CD11c Monoclonal Antibody, FITC Conjugated Clone BU15 853.381.010

Product Name	Anti-Human CD11c Monoclonal Antibody, FITC Conjugated Clone BU15
Catalog Number	853.381.010
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	FITC
Tested Applications	IHC, FC, IP
Immunogen	Dendritic cells of synovial fluid
Brand	Diaclone
Form/Appearance	The reagent is provided in phosphate buffered saline (PBS) containing 15 mM sodium azide and 0.2% (w/v) high-grade protease free Bovine Serum Albumin (BSA) as a stabilizing agent
Storage	Store at +2-8°C protected from light. DO NOT FREEZE.
Note	The reagent is designed for Flow Cytometry analysis of human blood cells using 20 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension.
Isotype	IgG1
Clonality	Monoclonal
Specificity	The antibody BU15 reacts with CD11c (alphaX, p150), a 150 kDa integrin expressed mainly on dendritic cells and tissue macrophages.