

Product Datasheet

Anti-Human CD11c Monoclonal Antibody, Unconjugated Clone BU15 853.383.020

Product Name	Anti-Human CD11c Monoclonal Antibody, Unconjugated Clone BU15
Catalog Number	853.383.020
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Purified
Tested Applications	IP, IHC, FC
Immunogen	Dendritic cells of synovial fluid
Brand	Diaclone
Form/Appearance	Phosphate buffered saline (PBS) with 15 mM sodium azide, approx. pH 7.4
Storage	Store at +2-8°C. DO NOT FREEZE.
Note	We recommend that each laboratory determine an optimum working titre for use in its particular application
Isotype	IgG1
Clonality	Monoclonal
Specificity	The antibody BU15 reacts with CD11c (alphaX, p150), a 150 kDa integrin expressed mainly on dendritic cells and tissue macrophages.