Medix Biochemica

Oy Medix Biochemica Ab

Headquarters Klovinpellontie 3 FI-02180 Espoo

Finland

Email: sales@medixbiochemica.com

Medix Biochemica USA Inc.

10850 Metro Court Maryland Heights, MO 63043 United States of America

Email: sales@medixbiochemica.com

Product Datasheet

Anti-Human CD35 Monoclonal Antibody, Unconjugated **Clone UJ11** 853.403.020

Product Name

Anti-Human CD35 Monoclonal Antibody, Unconjugated Clone UJ11

Catalog 853.403.020

Number

Description Mouse monoclonal

Species/Host Mouse

Conjugation Purified

Tested Applications

IHC, FC, IP

Brand Diaclone

Form/AppearaPBS with 15 mM sodium azide (Ph7.4)

Store at 4°C. For long-term storage aliquot and store at -20°C. Storage

Avoid freeze/thaw cycles.

Note We recommend that each laboratory determine an optimum

working titre for use in its particular application

Isotype lgG1,k

Clonality Monoclonal

Specificity The CD35 antigen, CR1, is a complement component receptor

with specificity for C3b, C4b and iC3b