

# Product Datasheet

## Anti-Human CD35 Monoclonal Antibody, Unconjugated Clone UJ11 853.403.020

Product Name	Anti-Human CD35 Monoclonal Antibody, Unconjugated Clone UJ11
Catalog Number	853.403.020
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Purified
Tested Applications	FC, IP, IHC
Brand	Diaclone
Form/Appearance	PBS with 15 mM sodium azide (Ph7.4)
Storage	Store at 4°C. For long-term storage aliquot and store at -20°C. Avoid freeze/thaw cycles.
Note	We recommend that each laboratory determine an optimum working titre for use in its particular application
Isotype	IgG1,k
Clonality	Monoclonal
Specificity	The CD35 antigen, CR1, is a complement component receptor with specificity for C3b, C4b and iC3b