

Product Datasheet

Semen Swab, Single Donor 991-04-S-SWB

Product Name	Semen Swab, Single Donor
Catalog Number	991-04-S-SWB
Description	Each sample is equivalent to 2 swabs per donor. Inquire for larger donor counts, swab counts, or custom criteria.
Tested Applications	QC Materials
Alternative Names	Seminal fluid, sperm, spermatoza, ejaculate, seed, germ, jism, male, male reproductive, reproductive, reproduction, forensic
Brand	Lee Biosolutions
Form/Appearance	Swabs
Storage	Ambient
Source	Single Human Donor
Shelf Life	Due to differing stabilities of various biomarkers in this complex biological material, it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time.
Preservatives	None added, inquire if needed