

Product Datasheet

Semen Swab - Normal Human Donor 991-04-S-SWB

Product Name	Semen Swab - Normal Human Donor
Catalog Number	991-04-S-SWB
Brand	Lee Biosolutions
Form/Appearance	Swab
Storage	-20 °C
Source	Single Human Donor
Shelf Life	Due to differing stabilities of various biomarkers in this complex biological material it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time.